

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Collagen I for tissue engineering - Bovine Skin, Unconjugated ROC-001-G70-MF7

Article Name	Collagen I for tissue engineering - Bovine Skin, Unconjugated
Biozol Catalog Number	ROC-001-G70-MF7
Supplier Catalog Number	001-G70-MF7
Alternative Catalog Number	ROC-001-G70-MF7
Manufacturer	Rockland Immunochemicals
Category	Proteine/Peptide
Application	SDS-PAGE
Conjugation	Unconjugated
Concentration	2.0 mg/mL by dry weight
Buffer	0.013M HCI
Form	Liquid (sterile filtered)
Application Notes	Bovine Collagen Type I solution has been tested in SDS-Page and can be utilized to mimic in vivo environment for 2D and 3D cell culture, tissue engineering research, and biochemistry.