



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## **TrueBlot Biotin Magnetic Beads ROC-B000-18-2**

Article Name	TrueBlot Biotin Magnetic Beads
Biozol Catalog Number	ROC-B000-18-2
Supplier Catalog Number	B000-18-2
Alternative Catalog Number	ROC-B000-18-2
Manufacturer	Rockland Immunochemicals
Category	Proteine/Peptide
Conjugation	Biotin
Concentration	5.0 mg/mL
Buffer	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Form	Suspension of Magnetic Beads
Application Dilute	IP: User Optimized, WB: User Optimized
Application Notes	TrueBlot Biotin magnetic beads are stored in phosphate buffered saline with detergent. TrueBlot Biotin magnetic beads can be used in immunoassays or molecular diagnostics when avidin or streptavidin is bound to proteins or oligonucleotides. Biotin magn