

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### **Bovine Serum Albumin 30% Solution, Unconjugated ROC-BSA-30-0010**

Article Name	Bovine Serum Albumin 30% Solution, Unconjugated
Biozol Catalog Number	ROC-BSA-30-0010
Supplier Catalog Number	BSA-30-0010
Alternative Catalog Number	ROC-BSA-30-0010
Manufacturer	Rockland Immunochemicals
Category	Zellen/Zellkultur
Conjugation	Unconjugated
Concentration	30% by dry weight
Buffer	0.85% (w/v) Sodium Chloride
Form	Liquid (sterile filtered)
Application Dilute	ELISA: User Optimized, FLISA: User Optimized, IHC: User Optimized, IF Microscopy: User Optimized, IP: User Optimized, WB: User Optimized
Application Notes	BOVINE SERUM ALBUMIN 30% Solution is suitable for use in protease sensitive assays such as RIA, EIA and nucleic acid hybridization, use as a stabilizing agent for proteins and enzymes, including dilute solutions of antibody, and use as a blocking agent t