

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Protein G Magnetic Beads IP/Co-IP Kit ROC-KBA-PG00-18

|                            |   |
|----------------------------|---|
| Article Name               | Protein G Magnetic Beads IP/Co-IP Kit   |
| Biozol Catalog Number      | ROC-KBA-PG00-18   |
| Supplier Catalog Number    | KBA-PG00-18   |
| Alternative Catalog Number | ROC-KBA-PG00-18   |
| Manufacturer               | Rockland Immunochemicals  |
| Category                   | Kits/Assays   |
| Application                | IP  |
| Concentration              | 5.0mg/mL by UV absorbance at 280 nm   |
| Form                       | Suspension of Magnetic Beads  |
| Application Dilute         | IP: User Optimized, WB: User Optimized  |
| Application Notes          | TrueBlot Protein G magnetic beads IP/Co-IP kit can be used for immunoprecipitation (IP) and Co-IP experiments as well as for antibody purification. For IP, target specific antibody is incubated with cell lysate. Protein G magnetic beads are then incubat |