

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Histone Targeted ChIP-seq Kit ROC-KOA0880

Article Name	Histone Targeted ChIP-seq Kit
Biozol Catalog Number	ROC-KOA0880
Supplier Catalog Number	KOA0880
Alternative Catalog Number	ROC-KOA0880
Manufacturer	Rockland Immunochemicals
Category	Kits/Assays
Application Dilute	ChIP: User Optimized
Application Notes	The different steps of the ChIP assay are cell fixation (cross-linking), chromatin shearing, immunoprecipitation, reverse cross-linking followed by DNA purification and analysis of the immunoprecipitated DNA. Suitable for 10 reactions.