



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## **Genomic DNA Magnetic Bead Kit ROC-KOA0912**

Article Name	Genomic DNA Magnetic Bead Kit
Biozol Catalog Number	ROC-KOA0912
Supplier Catalog Number	KOA0912
Alternative Catalog Number	ROC-KOA0912
Manufacturer	Rockland Immunochemicals
Category	Kits/Assays
Concentration	12.5 mg/mL
Form	Suspension of Magnetic Beads
Application Dilute	IP: User Optimized, WB: User Optimized
Application Notes	The Genomic DNA Magnetic Beads Kit can be used for extraction and isolation of genomic DNA from bacteria and mammalian cells and tissue samples.