



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## **Genomic DNA Magnetic Bead Kit ROC-KOA0912**

| Article Name               | Genomic DNA Magnetic Bead Kit                                                                                                                    |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------|
| Biozol Catalog Number      | ROC-KOA0912                                                                                                                                      |
| Supplier Catalog Number    | KOA0912                                                                                                                                          |
| Alternative Catalog Number | ROC-KOA0912                                                                                                                                      |
| Manufacturer               | Rockland Immunochemicals                                                                                                                         |
| Category                   | Kits/Assays                                                                                                                                      |
| Concentration              | 12.5 mg/mL                                                                                                                                       |
| Form                       | Suspension of Magnetic Beads                                                                                                                     |
| Application Dilute         | IP: User Optimized, WB: User Optimized                                                                                                           |
| Application Notes          | The Genomic DNA Magnetic Beads Kit can be used for extraction and isolation of genomic DNA from bacteria and mammalian cells and tissue samples. |