



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Collagenase Type 2, Unconjugated ROC-MB-119-0100

Article Name	Collagenase Type 2, Unconjugated
	conagenase Type 2, oneonjagatea
Biozol Catalog Number	ROC-MB-119-0100
Supplier Catalog Number	MB-119-0100
Alternative Catalog Number	ROC-MB-119-0100
Manufacturer	Rockland Immunochemicals
Category	Proteine/Peptide
Conjugation	Unconjugated
Concentration	125 CDU/mgdw by dry weight
NCBI	46085
UniProt	Q46085
Buffer	None
Source	Clostridium histolyticum
Form	Lyophilized
Application Notes	Collagenase Type 2 contain higher clostripain activity. This is suggested for bone, heart, liver, thyroid and salivary primary cell isolation. Collagenase is typically used at concentrations from 0.05 % to 0.5 % (w/v) in balanced salt solutions such as H