

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Protein Molecular Weight Marker, Unconjugated ROC-MB-201-1000

Article Name	Protein Molecular Weight Marker, Unconjugated
Biozol Catalog Number	ROC-MB-201-1000
Supplier Catalog Number	MB-201-1000
Alternative Catalog Number	ROC-MB-201-1000
Manufacturer	Rockland Immunochemicals
Category	Proteine/Peptide
Application	SDS-PAGE
Conjugation	Unconjugated
Concentration	1X
Buffer	2.5 mM Tris-HCl, pH 6.8, 2% (w/v) SDS, 250 mM NaCl, 1mM EDTA, 4% (v/v) beta-Mercaptoethanol, 35% (v/v) Glycerol and 0.05% (w/v) bromophenol blue
Form	Liquid (in 1x Loading Buffer)
Application Notes	Protein Molecular Weight Marker has been tested and is designed for use in accurate estimation of protein molecular weight in SDS-polyacrylamide gel electrophoresis and western blotting. Thaw the cocktail at room temperature. Gently invert vial several t