

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## Opal Prestained Protein Standard 10-180kDa, Unconjugated ROC-MB-209-0500

Article Name	Opal Prestained Protein Standard 10-180kDa, Unconjugated
Biozol Catalog Number	ROC-MB-209-0500
Supplier Catalog Number	MB-209-0500
Alternative Catalog Number	ROC-MB-209-0500
Manufacturer	Rockland Immunochemicals
Category	Proteine/Peptide
Application	SDS-PAGE, WB
Conjugation	Unconjugated
Concentration	1X
Buffer	See application note.
Form	Liquid (in 1x Loading Buffer)
Application Notes	Opal Prestained Protein Standard 10-180kDa is supplied in gel loading buffer and is ready to use. Do not heat, dilute, add reducing agent before loading. Suggested load volumes: 3 $\mu$ l or 5 $\mu$ l per loading for clear visualization during electrophoresis on 1