

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Protein G Peroxidase Conjugated, HRP ROC-PG00-03-0500

Article Name	Protein G Peroxidase Conjugated, HRP
Biozol Catalog Number	ROC-PG00-03-0500
Supplier Catalog Number	PG00-03-0500
Alternative Catalog Number	ROC-PG00-03-0500
Manufacturer	Rockland Immunochemicals
Category	Proteine/Peptide
Application	DOT, ELISA
Conjugation	HRP
Concentration	1.0 mg/mL by UV absorbance at 280 nm
NCBI	<a href="#">19909</a>
UniProt	<a href="#">P19909</a>
Buffer	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Form	Lyophilized
Application Dilute	ELISA: 1:20,000 - 1:200,000, IHC: 1:1,000 - 1:5,000, WB: 1:10,000 - 1:40,000
Application Notes	Protein G Peroxidase has been tested by ELISA and dot blot and is a useful reagent in Western Blotting and ELISA experiments. Protein G Peroxidase can be utilized as a pseudo-secondary detection reagent when used in conjunction with an IgG-based primary