



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Protein G Peroxidase Conjugated, HRP ROC-PG00-03-0500

| Article Name | Protein G Peroxidase Conjugated, HRP |
|----------------------------|---|
| Biozol Catalog Number | ROC-PG00-03-0500 |
| Supplier Catalog Number | PG00-03-0500 |
| Alternative Catalog Number | ROC-PG00-03-0500 |
| Manufacturer | Rockland Immunochemicals |
| Category | Proteine/Peptide |
| Application | DOT, ELISA |
| Conjugation | HRP |
| Concentration | 1.0 mg/mL by UV absorbance at 280 nm |
| NCBI | 19909 |
| UniProt | P19909 |
| Buffer | 0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2 |
| Form | Lyophilized |
| Application Dilute | ELISA: 1:20,000 - 1:200,000, IHC: 1:1,000 - 1:5,000, WB: 1:10,000 - 1:40,000 |
| Application Notes | Protein G Peroxidase has been tested by ELISA and dot blot and is a useful reagent in Western Blotting and ELISA experiments. Protein G Peroxidase can be utilized as a pseudo-secondary detection reagent when used in conjunction with an IgG-based primary |