



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

SEPHAROSE(TM) PROTEIN G, Sepharose ROC-PG50-00-0002

Article Name	SEPHAROSE(TM) PROTEIN G, Sepharose
Biozol Catalog Number	ROC-PG50-00-0002
Supplier Catalog Number	PG50-00-0002
Alternative Catalog Number	ROC-PG50-00-0002
Manufacturer	Rockland Immunochemicals
Category	Proteine/Peptide
Conjugation	Sepharose
NCBI	19909
UniProt	P19909
Form	Suspension of agarose beads
Application Dilute	IP: User Optimized
Application Notes	Protein G resin may also be used to pull down antibody:antigen complexes in immunoprecipitation experiments. This product is suitable for use at 20 μ L per immunoprecipitation reaction.