



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## Glycine max (SBA) (Colloidal Gold) USB-518373

Article Name	Glycine max (SBA) (Colloidal Gold)
Biozol Catalog Number	USB-518373
Supplier Catalog Number	518373
Alternative Catalog Number	USB-518373-1
Manufacturer	US Biological
Category	Molekularbiologie
Product Description	Colloidal gold can be used to visualize localization of specific saccharide moieties and examine the distribution of plasma membrane glycoproteins based on the binding pattern of Glycine max (Soybean) Lectin (SBA). Colloidal gold labeled lectins can
Purity	Highly Purified
Form	Supplied as a liquid in 10mM HEPES, 150mM NaCl, 0.1mM lons, 0.02% sodium azide. Labeled with Gold.