



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## Glycine max (SBA) (DyLight488) USB-518376

Article Name	Glycine max (SBA) (DyLight488)
Biozol Catalog Number	USB-518376
Supplier Catalog Number	518376
Alternative Catalog Number	USB-518376-1
Manufacturer	US Biological
Category	Molekularbiologie
Product Description	Dylight 488 can be used to visualize the binding pattern of Glycine max (Soybean) Lectin (SBA) in cellular imaging and flow cytometry. DyLight 488 is a green-fluorescent dye that works over a broad pH range (pH 4-9). It is brighter and more photostab
Form	Supplied as a liquid in 10mM HEPES, 150mM NaCl, 0.1mM lons, 0.02% sodium azide. Labeled with DyLight488.