



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## Glycine max (SBA) (FITC) USB-518378

Article Name	Glycine max (SBA) (FITC)
Biozol Catalog Number	USB-518378
Supplier Catalog Number	518378
Alternative Catalog Number	USB-518378-2,USB-518378-5
Manufacturer	US Biological
Category	Molekularbiologie
Product Description	Glycine max (Soybean) Lectin (SBA) is labeled with fluorescein isothiocyanate (FITC) and has an appropriate number of fluorochromes bound to provide the optimum staining characteristics for this lectin. FITC conjugates have been used in a variety of
Purity	Affinity Purified
Form	Supplied as a liquid in 10mM HEPES, 150mM NaCl, 0.1mM lons, 0.02% sodium azide. Labeled with Fluorescein isothiocyanate (FITC).