



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Lens culinaris (LCA/LCH) (Colloidal Gold) USB-518462

Article Name	Lens culinaris (LCA/LCH) (Colloidal Gold)
Biozol Catalog Number	USB-518462
Supplier Catalog Number	518462
Alternative Catalog Number	USB-518462-1
Manufacturer	US Biological
Category	Molekularbiologie
Product Description	Colloidal gold can be used to visualize localization of specific saccharide moieties and examine the distribution of plasma membrane glycoproteins based on the binding pattern of Lens culinaris (Lentil) Lectin (LCA/LCH). Colloidal gold labeled lectin
Purity	Highly Purified
Form	Supplied as a liquid in 10mM HEPES, 150mM NaCl, 0.1mM lons, 0.02% sodium azide. Labeled with Gold.