



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## Lotus tetragonolobus (LTL) (Cy3) USB-518481

Article Name	Lotus tetragonolobus (LTL) (Cy3)
Biozol Catalog Number	USB-518481
Supplier Catalog Number	518481
Alternative Catalog Number	USB-518481-1
Manufacturer	US Biological
Category	Molekularbiologie
Product Description	Cy3 can be used to viualize the binding pattern of Lotus tetragonolobus Lectin (LTL) in cellular imaging and flow cytometry. Cy3 is more photostable than many other flourophores and can be seen with TRITC filter sets. It is commonly combined with gre
Form	Supplied as a liquid in 10mM HEPES, 150mM NaCl, 0.1mM lons, 0.02% sodium azide. Labeled with Cy3.