



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Phaseolus vulgaris (PHA-E+L) (Cy3) USB-518637

Article Name	Phaseolus vulgaris (PHA-E+L) (Cy3)
Biozol Catalog Number	USB-518637
Supplier Catalog Number	518637
Alternative Catalog Number	USB-518637-1
Manufacturer	US Biological
Category	Molekularbiologie
Product Description	Cy3 can be used to visualize the binding pattern of Phaseolus vulgaris Lectin (PHA-E+L) in cellular imaging and flow cytometry. Cy3 is more photostable than many other fluorophores and can be seen with TRITC filter sets. It is commonly combined with
Form	Supplied as a liquid in 10mM HEPES, 150mM NaCl, 0.1mM lons, 0.02% sodium azide. Labeled with Cy3.