



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## Solanum tuberosum (STA/STL) (Cy3) USB-518757

Article Name	Solanum tuberosum (STA/STL) (Cy3)
Biozol Catalog Number	USB-518757
Supplier Catalog Number	518757
Alternative Catalog Number	USB-518757-1,USB-518757-2
Manufacturer	US Biological
Category	Molekularbiologie
Product Description	Cy3 can be used to visualize the binding pattern of Solanum tuberosum Lectin (STA/STL) in cellular imaging and flow cytometry. Cy3 is more photostable than many other fluorophores and can be seen with TRITC filter sets. It is commonly combined with g
Form	Supplied as a liquid in PBS, 0.02% sodium azide. Labeled with Cy3.