



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Histone, acetylated BioAssay(TM) Sample Kit USB-H5110-05

Article Name	Histone, acetylated BioAssay(TM) Sample Kit
Biozol Catalog Number	USB-H5110-05
Supplier Catalog Number	H5110-05
Alternative Catalog Number	USB-H5110-05-8X20
Manufacturer	US Biological
Category	Kits/Assays
Product Description	The Acetyl-Histone Antibody Sampler Kit provides a fast and economical means to evaluate the activation status of the acetylation sites on all of the histones. The kit contains enough primary and secondary antibodies to perform two mini-blot experime
Purity	Purified by Protein A and peptide affinity chromatography.
Form	Supplied as a liquid in 10mM sodium HEPES, pH 7.5, 150mM sodium chloride, 0.1mg/ml BSA, 50% glycerol.