



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Histone H3, phosphorylated (Ser10) Sandwich ELISA Kit, BioAssay(TM) USB-H5110-12N

Article Name	Histone H3, phosphorylated (Ser10) Sandwich ELISA Kit, BioAssay(TM)
Biozol Catalog Number	USB-H5110-12N
Supplier Catalog Number	H5110-12N
Alternative Catalog Number	USB-H5110-12N-1
Manufacturer	US Biological
Category	Kits/Assays
Product Description	Detects endogenous levels of Phospho-Histone H3 (Ser10). Using this Sandwich ELISA Kit Phospho- Histone H3 (Ser10) is detected when treated with Calyculin A in NIH 3T3 cells. However, the levels of Histone H3 remains unchanged, as shown by Western an