

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## In situ Hybridization Kit BioAssay(TM)(ISH), Blocking Solution, 100X USB-I7643-50H

Article Name	In situ Hybridization Kit BioAssay(TM)(ISH), Blocking Solution, 100X
Biozol Catalog Number	USB-I7643-50H
Supplier Catalog Number	17643-50H
Alternative Catalog Number	USB-I7643-50H-250
Manufacturer	US Biological
Category	Kits/Assays
Product Description	100X Blocking Solution for I7643-50: In situ Hybridization Kit (ISH), BioAssay(TM) Storage and Stability: Store at 4C Stable for 6 months after receipt For maximum recovery of product, centrifuge the original vial prior to removing the cap