

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## Kinase Substrate, phosphorylated (Ser/Thr) BioAssay(TM) Sample Kit USB-K0160

Article Name	Kinase Substrate, phosphorylated (Ser/Thr) BioAssay(TM) Sample Kit
Biozol Catalog Number	USB-K0160
Supplier Catalog Number	K0160
Alternative Catalog Number	USB-K0160-6X20
Manufacturer	US Biological
Category	Kits/Assays
Product Description	Phospho-(Ser/Thr) kinases and phosphatases play critical roles in a wide range of biological processes. Each phospho-(Ser/Thr) kinase phosphorylates serine or threonine within a specific motif. Akt phosphorylates substrates at a serine or threonine o
Purity	Purified by Protein A and peptide affinity chromatography.
Form	Supplied as a liquid in 10mM sodium HEPES, pH 7.5, 150mM sodium chloride, 0.1mg/ml BSA, 50% glycerol.