



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

MAPK Family BioAssay(TM) Sample Kit USB-M2352-40

Article Name	MAPK Family BioAssay(TM) Sample Kit
Biozol Catalog Number	USB-M2352-40
Supplier Catalog Number	M2352-40
Alternative Catalog Number	USB-M2352-40-1
Manufacturer	US Biological
Category	Kits/Assays
Application	ELISA, WB
Product Description	P44/42 MAPK (Erk1 and Erk2), SAPK and p38 MAPK function in protein kinase cascades that play a critical role in the regulation of cell growth, and the differentiation and control of cellular responses to cytokines and stress. p44/42 MAPK (Erk1 and Er
Purity	Purified by Protein A and peptide affinity chromatography.
Form	Supplied as a liquid in 10mM HEPES, pH 7.5, 150mM sodium chloride, 0.1mg/ml BSA, 50% glycerol. No preservative added.