

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

p-Nitrophenyl Phosphate (pNPP) Tyrosine Assay Buffer USB-N2600-30

Article Name	p-Nitrophenyl Phosphate (pNPP) Tyrosine Assay Buffer
Biozol Catalog Number	USB-N2600-30
Supplier Catalog Number	N2600-30
Alternative Catalog Number	USB-N2600-30-20
Manufacturer	US Biological
Category	Kits/Assays
Product Description	Applications: For use in the Tyrosine Assay Kit. Inclusion of dithiothreitol is critical for optimal PTP activity. Storage and Stability: May be stored at 4C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at 4C. Aliquots a
Purity	Sterilized through a 0.2um filter and packaged aseptically
Form	Supplied as a liquid in 25mM HEPES, pH 7.2, 50mM sodium chloride, 5mM dithiothreitol, 2.5mM EDTA.