



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Nitric Oxide Synthase Kit, BioAssay(TM) (NOS) USB-N5350-33

Article Name	Nitric Oxide Synthase Kit, BioAssay(TM) (NOS)
Biozol Catalog Number	USB-N5350-33
Supplier Catalog Number	N5350-33
Alternative Catalog Number	USB-N5350-33-96
Manufacturer	US Biological
Category	Kits/Assays
Product Description	The traditional method for measuring nitric oxide synthase (NOS) activity is performed by radiochemical assay that measures the conversion of L-[3H]arginine to L-[3H]citrulline. This method is expensive and requires regulation of radioactive material