



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Presporulation Medium (Powder) USB-P8900

Article Name	Presporulation Medium (Powder)
Biozol Catalog Number	USB-P8900
Supplier Catalog Number	P8900
Alternative Catalog Number	USB-P8900-500,USB-P8900-2.5,USB-P8900-10
Manufacturer	US Biological
Category	Biochemikalien
Product Description	Presporulation Medium is used on strains of yeast cells that do not sporulate efficiently, such as those that carry mutations that cause the cells to grow slowly or ferment poorly, are grown for 1 to 2 days before the cells are transferred to Sporula