



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## 2-Thiobarbituric Acid Reactive Substances (TBARS) Colorimetric BioAssay(TM) Kit USB-T0700-89

Article Name	2-Thiobarbituric Acid Reactive Substances (TBARS) Colorimetric BioAssay(TM) Kit
Biozol Catalog Number	USB-T0700-89
Supplier Catalog Number	T0700-89
Alternative Catalog Number	USB-T0700-89-40
Manufacturer	US Biological
Category	Kits/Assays
Product Description	The TBARS assay was developed for the measurement of MDA production via lipid peroxidation and is a well-established and widely used method for assessment of oxidative damage. However, as typically performed, the traditional method for TBARS analysis