

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Translational Control BioAssay(TM) Sample Kit USB-T8263

Article Name	Translational Control BioAssay(TM) Sample Kit
Biozol Catalog Number	USB-T8263
Supplier Catalog Number	T8263
Alternative Catalog Number	USB-T8263-6X20
Manufacturer	US Biological
Category	Kits/Assays
Product Description	The synthesis of new protein is a highly regulated process that allows rapid cellular responses to diverse stimuli in the absence of transcription. Initiation of protein synthesis begins after the separation of the ribosome into its 40S and 60S subun
Purity	Purified by Protein A and peptide affinity chromatography.
Form	Supplied as a liquid in 10mM sodium HEPES, pH 7.5, 150mM sodium chloride, 0.1mg/ml BSA, 50% glycerol.