



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Recombinant Histone H4 (Xenopus), E. coli MBS388037-5X0.05

Article Name	Recombinant Histone H4 (Xenopus), E. coli
Biozol Catalog Number	MBS388037-5X0.05
Supplier Catalog Number	MBS388037-5x0.05
Alternative Catalog Number	MBS388037-5X0.05
Manufacturer	MyBiosource
Host	E. coli
Category	Proteine/Peptide
Purity	The recombinant histone is >98% pure by SDS-PAGE.
Form	Recombinant Histone H4 (Xenopus) is supplied as lyophilized powder. Recombinant histones can be resuspended in water or any suitable buffer. We recommend a starting concentration of 1mg/ml. To fully solubilize the histone we suggest resuspension in the b