



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Recombinant Histone H3K23ac (EPL), E. coli MBS388043-5X0.025

| Article Name | Recombinant Histone H3K23ac (EPL), E. coli |
|----------------------------|--|
| Biozol Catalog Number | MBS388043-5X0.025 |
| Supplier Catalog Number | MBS388043-5x0.025 |
| Alternative Catalog Number | MBS388043-5X0.025 |
| Manufacturer | MyBiosource |
| Host | E. coli |
| Category | Proteine/Peptide |
| Purity | The recombinant histone is >98% pure by SDS-PAGE and confirmed by high-resolution ESI-TOF mass spectrometry |
| Form | 25ug supplied as lyophilized powder. Recombinant histones can be resuspended in water or any suitable buffer. We recommend a starting concentration of 1mg/ml. To fully solubilize the histone we suggest resuspension in the buffer of choice at room tempera |