



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## PurKine(TM) GST-Tag Glutathione Resin ASC-BMR2010

Article Name	PurKine(TM) GST-Tag Glutathione Resin
Biozol Catalog Number	ASC-BMR2010
Supplier Catalog Number	BMR2010
Alternative Catalog Number	ASC-BMR2010-5,ASC-BMR2010-25
Manufacturer	Abbkine Scientific
Category	Sonstiges
Product Description	PurKine(TM) GST Tagged Protein Purification Glutathione Resin
Tag	GST Tag
Buffer	Liquid in PBS buffer containing 20% ethanol.
Form	Liquid solution with 50% slurry in 20% ethanol.
Application Notes	PurKine(TM) GST-Tag Glutathione Resin effectively purifies high levels of overexpressed GST-tagged fusion proteins in gravity column procedures at a variety of scales. The Resin consists of 90µm beads of cross-linked 4% agarose, to which reduced glutathione