

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

PurKine(TM) Strep II-Tag Strep-Tactin Resin 4FF ASC-BMR2040

Article Name	PurKine(TM) Strep II-Tag Strep-Tactin Resin 4FF
Biozol Catalog Number	ASC-BMR2040
Supplier Catalog Number	BMR2040
Alternative Catalog Number	ASC-BMR2040-2,ASC-BMR2040-10
Manufacturer	Abbkine Scientific
Category	Sonstiges
Product Description	PurKine(TM) Strep II-Tag Strep-Tactin Resin 4FF
Tag	Strep II-Tag
Buffer	Liquid in PBS buffer containing 20% ethanol.
Form	Liquid solution with 50% slurry in 20% ethanol.
Application Notes	PurKine(TM) Strep II-Tag Strep-Tactin Resin 4FF effectively purifies high levels of overexpressed Strep II-tagged fusion proteins in gravity column procedures at a variety of scales. The Resin consists of 90µm beads of cross-linked 4% agarose, to which Stre