

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

SuperKine(TM) Maximum Sensitivity Cell Counting Kit-8 (CCK-8) ASC-BMU106-EN

Article Name	SuperKine(TM) Maximum Sensitivity Cell Counting Kit-8 (CCK-8)
Biozol Catalog Number	ASC-BMU106-EN
Supplier Catalog Number	BMU106-EN
Alternative Catalog Number	ASC-BMU106-EN-200,ASC-BMU106-EN-1000,ASC-BMU106-EN-10000
Manufacturer	Abbkine Scientific
Category	Zellen/Zellkultur
Product Description	SuperKine(TM) Maximum Sensitivity Cell Counting Kit-8 (CCK-8) allows very convenient assays by utilizing highly water-soluble tetrazolium salt-WST-8. WST-8 is reduced by dehydrogenases in cells to give an orange colored product (Formazan), which is s...
Tag	Cell proliferation/toxicity was detected by CCK-8
Form	Liquid solution
Application Notes	Wst-8 was an upgrade of MTT and was reduced by mitochondrial dehydrogenase in the presence of electron-coupled reagents to form a highly water-soluble orange dirty product (Formazan).