



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Human Tapasin (TAPBP) ELISA Kit ASC-KTE60473

| Article Name | Human Tapasin (TAPBP) ELISA Kit |
|----------------------------|--|
| Biozol Catalog Number | ASC-KTE60473 |
| Supplier Catalog Number | KTE60473 |
| Alternative Catalog Number | ASC-KTE60473-48,ASC-KTE60473-96 |
| Manufacturer | Abbkine Scientific |
| Category | Kits/Assays |
| Application | ELISA |
| Species Reactivity | Human |
| Product Description | This Human Tapasin (TAPBP) ELISA Kit employs a two-site sandwich ELISA to quantitate TAPBP |
| Range | Please inquire |
| Sensitivity | Please inquire |
| Tag | TAPBP |
| NCBI | 6892 |
| UniProt | O15533 |
| Samples | Cell culture supernatantsSerumPlasmaOther biological fluids |
| Application Notes | This Human Tapasin (TAPBP) ELISA Kit employs a two-site sandwich ELISA to quantitate TAPBP in samples. An antibody specific for TAPBP has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and anyTAPBP present is bound b |