

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## Human Brain natriuretic peptide (BNP) ELISA Kit ASC-KTE63076

| Article Name               | Human Brain natriuretic peptide (BNP) ELISA Kit  |
|----------------------------|--|
| Biozol Catalog Number      | ASC-KTE63076   |
| Supplier Catalog Number    | KTE63076   |
| Alternative Catalog Number | ASC-KTE63076-96,ASC-KTE63076-48  |
| Manufacturer               | Abbkine Scientific   |
| Category                   | Kits/Assays  |
| Application                | ELISA  |
| Species Reactivity         | Human  |
| Product Description        | Null employs a two-site sandwich ELISA to quantitate BNP in samples  |
| Range                      | 15 pg/mL-240 pg/mL   |
| Sensitivity                | 1.0 pg/mL  |
| Tag                        | BNP  |
| Samples                    | Cell culture supernatantsSerumPlasmaOther biological fluids  |
| Application Notes          | Human Brain natriuretic peptide (BNP) ELISA Kit employs a two-site sandwich ELISA to quantitate BNP in samples. An antibody specific for BNP has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any BNP present is b |