

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### IPKine(TM) Anti-GST Magnetic IP Kit, Magnetic Beads ASC-KTI2034

|                            |  |
|----------------------------|--|
| Article Name               | IPKine(TM) Anti-GST Magnetic IP Kit, Magnetic Beads  |
| Biozol Catalog Number      | ASC-KTI2034  |
| Supplier Catalog Number    | KTI2034  |
| Alternative Catalog Number | ASC-KTI2034-20,ASC-KTI2034-100   |
| Manufacturer               | Abbkine Scientific   |
| Category                   | Antikörper   |
| Application                | IP   |
| Species Reactivity         | Bacteria, Mammal   |
| Immunogen                  | Synthetic Peptide  |
| Conjugation                | Magnetic Beads   |
| Product Description        | The kit contains high-quality Anti-GST magnetic beads and the necessary reagents for immunoprecipitation that have been optimized and validated, making immunoprecipitation (IP) experiments simpler, more convenient and more efficient....                   |
| Tag                        | GST  |
| Application Notes          | Anti-Glutathione S-transferase-Tag (GST-Tag) Magnetic Beads are prepared by covalently coupling Anti-GST-Tag Mouse Monoclonal Antibody to crosslinked Magnetic Beads. The Optimized Anti-GST Magnetic Beads have more efficient antigen binding capacity. This |