

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### **LinKine(TM) AbFluor(TM) 488 Labeling Kit (Optimized for samples with molecular weight of 6 KD to 20 KD), AbFluor 488 ASC-KTL0521**

Article Name	LinKine(TM) AbFluor(TM) 488 Labeling Kit (Optimized for samples with molecular weight of 6 KD to 20 KD), AbFluor 488
Biozol Catalog Number	ASC-KTL0521
Supplier Catalog Number	KTL0521
Alternative Catalog Number	ASC-KTL0521-1
Manufacturer	Abbkine Scientific
Category	Kits/Assays
Conjugation	AbFluor 488
Product Description	LinKine(TM) AbFluor(TM) 488 Labeling Kit (Optimized for samples with molecular weight of 6 KD to 20 KD) is designed for preparing AbFluor(TM) 488 conjugates directly from proteins, peptides, and other ligands that contain a free amino group....
Tag	AbFluor(TM) 488
Application Notes	LinKine(TM) AbFluor(TM) 488 Labeling Kit (Optimized for samples with molecular weight of 6 KD to 20 KD) is designed for preparing AbFluor(TM) 488 conjugates directly from proteins, peptides, and other ligands that contain a free amino group.