

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### **LinKine(TM) AbFluor(TM) 555 Labeling Kit (Optimized for samples with molecular weight of 6 KD to 20 KD), AbFluor 555 ASC-KTL0531**

Article Name	LinKine(TM) AbFluor(TM) 555 Labeling Kit (Optimized for samples with molecular weight of 6 KD to 20 KD), AbFluor 555
Biozol Catalog Number	ASC-KTL0531
Supplier Catalog Number	KTL0531
Alternative Catalog Number	ASC-KTL0531-1
Manufacturer	Abbkine Scientific
Category	Kits/Assays
Conjugation	AbFluor 555
Product Description	LinKine(TM) AbFluor(TM) 555 Labeling Kit (Optimized for samples with molecular weight of 6 KD to 20 KD) is designed for preparing AbFluor(TM) 555 conjugates directly from proteins, peptides, and other ligands that contain a free amino group....
Tag	AbFluor(TM) 555
Application Notes	LinKine(TM) AbFluor(TM) 555 Labeling Kit (Optimized for samples with molecular weight of 6 KD to 20 KD) is designed for preparing AbFluor(TM) 555 conjugates directly from proteins, peptides, and other ligands that contain a free amino group.