

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### **LinKine(TM) AbFluor(TM) 647 Labeling Kit (Optimized for samples with molecular weight of 6 KD to 20 KD), AbFluor 647 ASC-KTL0561**

Article Name	LinKine(TM) AbFluor(TM) 647 Labeling Kit (Optimized for samples with molecular weight of 6 KD to 20 KD), AbFluor 647
Biozol Catalog Number	ASC-KTL0561
Supplier Catalog Number	KTL0561
Alternative Catalog Number	ASC-KTL0561-1
Manufacturer	Abbkine Scientific
Category	Kits/Assays
Conjugation	AbFluor 647
Product Description	LinKine(TM) AbFluor(TM) 647 Labeling Kit (Optimized for samples with molecular weight of 6 KD to 20 KD) is designed for preparing AbFluor(TM) 647 conjugates directly from proteins, peptides, and other ligands that contain a free amino group....
Tag	AbFluor(TM) 647
Application Notes	LinKine(TM) AbFluor(TM) 647 Labeling Kit (Optimized for samples with molecular weight of 6 KD to 20 KD) is designed for preparing AbFluor(TM) 647 conjugates directly from proteins, peptides, and other ligands that contain a free amino group.