

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

LinKine(TM) AbFluor(TM) 680 Labeling Kit (Optimized for samples with molecular weight of 6 KD to 20 KD), AbFluor 680 ASC-KTL0581

Article Name	LinKine(TM) AbFluor(TM) 680 Labeling Kit (Optimized for samples with molecular weight of 6 KD to 20 KD), AbFluor 680
Biozol Catalog Number	ASC-KTL0581
Supplier Catalog Number	KTL0581
Alternative Catalog Number	ASC-KTL0581-1
Manufacturer	Abbkine Scientific
Category	Kits/Assays
Conjugation	AbFluor 680
Product Description	LinKine(TM) AbFluor(TM) 680 Labeling Kit (Optimized for samples with molecular weight of 6 KD to 20 KD) is designed for preparing AbFluor(TM) 680 conjugates directly from proteins, peptides, and other ligands that contain a free amino group....
Tag	AbFluor(TM) 680
Application Notes	LinKine(TM) AbFluor(TM) 680 Labeling Kit (Optimized for samples with molecular weight of 6 KD to 20 KD) is designed for preparing AbFluor(TM) 680 conjugates directly from proteins, peptides, and other ligands that contain a free amino group.