

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Anti-Galectin 3/LGALS3 Antibody Picoband(TM), Rabbit, Polyclonal BOB-PB9081

Article Name	Anti-Galectin 3/LGALS3 Antibody Picoband(TM), Rabbit, Polyclonal
Biozol Catalog Number	BOB-PB9081
Supplier Catalog Number	PB9081
Alternative Catalog Number	BOB-PB9081-10UG,BOB-PB9081-100UG
Manufacturer	Boster Bio
Host	Rabbit
Category	Antikörper
Application	FC, ICC, IF, IHC, WB
Species Reactivity	Human
Immunogen	E.coli-derived human Galectin 3 recombinant protein (Position: K139-I250). Human Galectin 3 shares 88% and 84% amino acid (aa) sequences identity with mouse and rat Galectin 3, respectively.
Product Description	Boster Bio Anti-Galectin 3/LGALS3 Antibody Picoband(TM) catalog PB9081. Tested in Flow Cytometry, IF, IHC, ICC, WB applications. This antibody reacts with Human
Clonality	Polyclonal
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 $\mu g/ml$.
Molecular Weight	26152 MW
UniProt	P17931

Buffer	Each vial contains 4mg Trehalose, 0.9mg NaCl, 0.2mg Na2HPO4, 0.05mg NaN3.
Purity	Immunogen affinity purified.
Form	Lyophilized
Application Dilute	Western blot, 0.1 - $0.5\mu g/ml$, Human Immunohistochemistry (Paraffinembedded Section), 0.5 - $1\mu g/ml$, Human, By Heat Immunocytochemistry/Immunofluorescence, $5\mu g/ml$, Human Flow Cytometry, 1 - $3\mu g/1x10$ 6 cells, Human
Application Notes	WB: The detection limit for Galectin 3 is approximately 0.25ng/lane under reducing conditions. Tested Species: In-house tested species with positive results. By Heat: Boiling the paraffin sections in 10mM citrate buffer, pH6.0, for 20mins is required for