

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Anti-VIP peptides VIP Antibody, Rabbit, Polyclonal BOB-RP1108

Article Name	Anti-VIP peptides VIP Antibody, Rabbit, Polyclonal
Biozol Catalog Number	BOB-RP1108
Supplier Catalog Number	RP1108
Alternative Catalog Number	BOB-RP1108-10UG,BOB-RP1108-100UG
Manufacturer	Boster Bio
Host	Rabbit
Category	Antikörper
Application	IF, IHC
Species Reactivity	Human, Mouse, Rat
Immunogen	A synthetic peptide corresponding to a sequence in the middle region of human VIP, different from the related mouse and rat sequences by four amino acids.
Product Description	Boster Bio Anti-VIP peptides VIP Antibody catalog RP1108. Tested in IF, IHC applications. This antibody reacts with Human, Mouse, Rat
Clonality	Polyclonal
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 $\mu g/ml$.
Molecular Weight	19169 MW
UniProt	P01282

Buffer	Each vial contains 5mg BSA, 0.9mg NaCl, 0.2mg Na2HPO4, 0.05mg NaN3.
Purity	Immunogen affinity purified.
Form	Lyophilized
Application Dilute	Immunohistochemistry (Paraffin-embedded Section), 0.5-1μg/ml, Human, Mouse, Rat, By Heat Immunofluorescence, 5 μg/ml, Rat
Application Notes	Tested Species: In-house tested species with positive results. By Heat: Boiling the paraffin sections in 10mM citrate buffer, pH6.0, for 20mins is required for the staining of formalin/paraffin sections. Other applications have not been tested. Optimal d