

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## E. coli DNA Topoisomerase II (Gyrase) Assay Kit Plus-100 (enzyme included) PFD-DSA100KE

| Article Name               | E. coli DNA Topoisomerase II (Gyrase) Assay Kit Plus-100 (enzyme included)   |
|----------------------------|--|
| Biozol Catalog Number      | PFD-DSA100KE   |
| Supplier Catalog Number    | DSA100KE   |
| Alternative Catalog Number | PFD-DSA100KE   |
| Manufacturer               | ProFoldin  |
| Category                   | Molekularbiologie  |
| Product Description        | Antibacterial Drug Target   DNA topoisomerases such as bacterial<br>topoisomerase II (gyrase) convert relaxed circular DNA into<br>supercoiled DNA (DNA supercoiling reaction). The DNA<br>Topoisomerase II (Gyrase) Assay Kit is based on the principle that<br>the s |