



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## **SNAP tag Nanoselector Agarose, Camelus HBO-HA710165**

| Article Name               | SNAP tag Nanoselector Agarose, Camelus   |
|----------------------------|--|
| Biozol Catalog Number      | HBO-HA710165   |
| Supplier Catalog Number    | HA710165   |
| Alternative Catalog Number | HBO-HA710165   |
| Manufacturer               | HUABIO   |
| Host                       | Camelus  |
| Category                   | Antikörper   |
| Application                | ChIP, IP, MS, Purification   |
| Product Description        | SNAP tag Nanoselector Agarose have been specifically designed to bind SNAP tag-fusion proteins. SNAP tag Nanoselector Agarose is based on small high-affinity recombinant alpaca antibody fragments covalently coupled to the surface of Agarose. SNAP tag |
| Buffer                     | 20% EtOH   |
| Form                       | Agarose slurry   |
| Application Dilute         | Immunoprecipitation (IP)/Co-IP ,Mass spectrometry (MS) ,Enzyme activity measurements   |