



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## RFP Nanoselector Magnetic beads, Camelus HBO-HA710168

Article Name	RFP Nanoselector Magnetic beads, Camelus
Biozol Catalog Number	HBO-HA710168
Supplier Catalog Number	HA710168
Alternative Catalog Number	HBO-HA710168
Manufacturer	HUABIO
Host	Camelus
Category	Antikörper
Application	ChIP, IP, MS, Purification
Product Description	RFP Nanoselector Magnetic beads have been specifically designed to bind RFP-fusion proteins. RFP Nanoselector Magnetic beads is based on small high-affinity recombinant alpaca antibody fragments covalently coupled to the surface of Magnetic beads. RF
Buffer	50 % slurry in PBS containing 20 % Ethanol,
Application Dilute	Immunoprecipitation (IP)/Co-IP Mass spectrometry (MS) Enzyme activity measurements