

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

BriClone SF Hybridoma Cloning Medium **QBS-BR-002**

Article Name	BriClone SF Hybridoma Cloning Medium
Biozol Catalog Number	QBS-BR-002
Supplier Catalog Number	BR-002
Alternative Catalog Number	QBS-BR-002-100
Manufacturer	QED Bioscience
Category	Zellen/Zellkultur
Application	Cell Culture
Product Description	BriClone SF is a Serum Free additive for the cloning medium used in the post-fusion stages of hybridoma production and for improving the efficiency of hybridoma cells cloning....
Concentration	20X Solution
Pubmed	7873262
Application Notes	Thaw and add to the hybridoma cloning medium as a 5% v/v supplement. BriClone SF can be used for: Hybridoma growth post-fusion (refer to your own protocol) Hybridoma Cloning (refer to your own protocol)