

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

BriClone Hybridoma Cloning Medium QBS-BRI10000

Article Name	BriClone Hybridoma Cloning Medium
Biozol Catalog Number	QBS-BRI10000
Supplier Catalog Number	BRI10000
Alternative Catalog Number	QBS-BRI10000-100
Manufacturer	QED Bioscience
Category	Zellen/Zellkultur
Application	Cell Culture
Product Description	BriClone is an additive for the cloning medium used in the post-fusion stages of hybridoma production and for improving the efficiency of hybridoma cells cloning....
Concentration	20X Solution
Pubmed	7873262
Application Notes	Thaw and add to the hybridoma cloning medium as a 5% v/v supplement. BriClone SF can be used for: Hybridoma growth post-fusion (refer to your own protocol) Hybridoma Cloning (refer to your own protocol)