



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## Rapid 26S Proteasome Purification Kit-S UBP-J4310

Article Name	Rapid 26S Proteasome Purification Kit-S
Biozol Catalog Number	UBP-J4310
Supplier Catalog Number	J4310
Alternative Catalog Number	UBP-J4310-1
Manufacturer	UBPBio, LLC
Category	Sonstiges
Species Reactivity	Human
Product Description	The Rapid 26S Proteasome Purification Kit is designed to rapidly purify the endogenous 26S proteasome from cell or tissue extracts within 8 hours according to a published method (reference 1) with modifications. This approach uses the N-terminal ubiq
Buffer	20 mM Tris, 150 mM NaCl, 2 mM BME, 10% Glycerol
Source	E. coli recombinant